

DAPI staining for cell cycle analysis

Binds preferentially A-T base regions in DNA. Process between 1 - 10 million of cells

- 1) Spin the cells at 500g, 5 min, 4°C.
- 2) Wash twice with 5mL cold PBS 1X, spin 5 min, 500 g, at 4°C.
- 3) Resuspend carefully the cell pellet in 500µL of PBS 1X (make sure to obtain a single-cell suspension).
- 4) Add the cell suspension to a 15mL Falcon tube containing 4.5 mL of ice-cold 70 % EtOH while vortexing slightly.
- 5) Incubate at 4°C for at least 2 hours. Cells can be stored in fixative at 4°C for several days, for longer storage keep them at -20.
- 6) Day of the staining: centrifuge the ethanol-suspended cells 5min at 1000g. Remove the supernatant carefully.
- 7) Resuspend the cell pellet in 5ml PBS 1X, wait 15 min at RT. During this time, count the cells again and distribute in tubes the same cell number (1-2 million cells)
- 8) Spin the cells 5min at 1000g.
- 9) Resuspend the cell pellet in 300µL of DAPI/TritonX-100.
- 10) Incubate 30 minutes at RT protected from light.
- 11) Transfer the sample to appropriate tubes for acquisition at the cytometer. Read DAPI signal in parameter: 405-450/50, lineal mode. Acquire your samples in Low.

DAPI/Triton X-100 Solution:

add 10µL of 1mg/ml DAPI to 10ml of 0.1% (v/v) TritonX-100 (prepared in PBS).

Important tips:

- ❖ Keep the same ratio DAPI:Number of cells
- ❖ Prepare always an asynchronous culture (untreated) to fix the positions of the G1 and G2 peaks.